


PDS No. 6641xx/ 664160-TRI	<b>PRODUCT DATA SHEET</b>			Page 1 of 1
Revision 08	Petri Dish, 100 x 20 mm			
	Greiner Item-No. 6641xx / 664160-TRI			
Valid for Item-No.:	<b>664102</b>	<b>664160</b> (sterile)	<b>664161</b> (sterile)	<b>664160-TRI</b> (sterile)

1.	Description / Specification	
1.1	Description	Petri Dish with vents, heavy design 664102: standard 664160: CELLSTAR®, physical surface treatment, sterile 664161: sterile 664160-TRI: CELLSTAR®, physical surface treatment, sterile, triple packed
1.2	Dimensions	See Customer Drawing Weight: lid: 8.1 – 8.6 g dish: 10.3 – 10.8 g
1.3	Volume	Max. volume: 100 ml Working volume: 16 - 17 ml 664160, -160-TRI: growth area: 58 cm <sup>2</sup>
1.4	Material / Resin	Dish: PS (Polystyrene), free of heavy metal Lid: PS (Polystyrene), free of heavy metal
1.5	Colour	Dish: clear Lid: clear
1.6	Sterilization	664102: no 664160, -161: SAL 10 <sup>-3</sup> 664160-TRI: SAL 10 <sup>-6</sup>
1.7	Quality Control	Raw Material-Control: physical testing Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	664102, -161: general laboratory product for bacteriology to be used by qualified personnel in a laboratory environment. 664160, -160-TRI: general laboratory product for cell culture to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- For single use only - Elevations for ventilation of culture

2.	Features	
2.1	Basic features	664102, -161: - 664160, -160-TRI: free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	664102: n/a 664160, -160-TRI: 4 years 664161: 5 years
2.7	Other Information	-

3.	Packaging	664102, -160, -161	664160-TRI
3.1	Pieces / Bag	15	15 (triple packed)
3.2	Pieces / Box	360	180
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	
3.4	Other Information	664160, -160-TRI: Certificate of Quality to download	

4.	Other Information
4.1	664160 / 664160-TRI: Research use only. Not for diagnostics.

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 07	Date 5 September 2022	Date 14 September 2022	Date 15 September 2022	
Date 26.08.2015	Name S. Kaelberer	Name Dr. R. Daum	Name A. Illig	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.