


PDS No. 6512xx	<b>PRODUCT DATA SHEET</b>				Page 1 of 1
Revision 06	96 Well Microplate, PP, V-bottom, Chimney Well				
	Greiner Item-No. 6512xx				
Valid for Item-No.:	<b>651201</b>	<b>651203</b>	<b>651204</b>	<b>651205</b>	<b>651206</b>
	<b>651209</b>	<b>651261</b>			

1.	Description / Specification	
1.1	Description	PP Microplate, 96 well, solid V-bottom, chimney well, alphanumeric well coding 651201, -203, -204, -205, -206 -209: non-sterile 651261: sterile
1.2	Dimensions	See Customer Drawing
1.3	Volume per well	Total volume: 340 µl (mathematically calculated) Working volume: 50 - 335 µl
1.4	Material / Resin	PP (Polypropylene), free of heavy metal
1.5	Colour	651201, -261: translucent 651203: red 651204: blue 651205: green 651206: yellow 651209: black
1.6	Sterilization	651201, -203, -204, -205, -206 -209: no 651261: SAL 10 <sup>-3</sup>
1.7	Quality Control	<u>Raw Material-Control</u> : physical testing <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for the processing and storage of samples to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- For single use only - Sealable with adhesive films and heat sealer or with CapMats

2.	Features	
2.1	Basic features	All Item-No.: free of detectable DNase/RNase, human DNA and pyrogens. 651261: contents non-cytotoxic
2.2	Temperature range	For application: -196°C to +121°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	4.800 x g: Swinging-bucket rotor
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	651201, -203, -204, -205, -206 -209: n/a 651261: 5 years
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	10
3.2	Pieces / Box	100
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality to download

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 05	Date 17 April 2024	Date 17 April 2024	Date 17 April 2024	
Date 21.01.2022	Name S. Kaelberer	Name P. Wachter	Name Dr. C.-K. Chai	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.