PDS No. 50x870	PRODUCT DATA SHEET			Page 1 of 1	
Revision 06	Glass Cover Slides, siliconized			0	
	Greiner Item-No. 50x870			greiner	
Valid for Item-No.:	501870	503870			

1.	Description / Specification		
1.1	Description	Glass cover slides, round, surface treated	
1.2	Dimensions	See Customer Drawing	
1.3	Volume	-	
1.4	Material / Resin	Cover: borosilicate glass	
		<u>Treatment</u> : siliconized	
1.5	Colour	Clear	
1.6	Sterilisation	No	
1.7	Quality Control	Product-Control: testing of attributive and variable characteristics in accordance with the valid specification	
1.8	Intended Use	General laboratory product for covering or capping to be used by qualified personnel in a laboratory environment.	
1.9	Other Information	For single use only	

2.	Features	
2.1	Basic features	-
2.2	Temperature range	N/A
2.3	Autoclavability	N/A
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	High resistant
2.6	Shelf life	N/A
2.7	Other Information	Production of glass cover slides in accordance with DIN ISO 8255

3.	Packaging	501870	503870
3.1	Pieces / Bag	50	100
3.2	Pieces / Box	1.000	1.000
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-	-

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this
Revision	Date	Date	Date	document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This
05	31 January 2023	2 February 2023		document may not be reproduced for any
Date	Name	Name	Name	reason without written permission from Greiner Bio-One GmbH. All rights of design, invention,
20.04.2015	S. Kaelberer	P. Wachter	A. Illig	and copyright are reserved.