


PDS No. 6091xx	<b>PRODUCT DATA SHEET</b>				Page 1 of 1
Revision 05	96 Well CrystalQuick / CrystalQuick Plus, PS				
	Greiner Item-No. 6091xx				
Valid for Item-No.:	<b>609101</b>	<b>609120</b>	<b>609130</b>	<b>609171</b>	<b>609180</b>

1.	Description / Specification	
1.1	Description	<p>Crystallisationplates, 96 well, with additional alphanumeric well coding per reservoir. For sitting drop applications with different well profiles.</p> <p>609101: CrystalQuick: square well design, flat bottom, 3 wells per reservoir</p> <p>609120: CrystalQuick: round well design, concave bottom, 3 wells per reservoir</p> <p>609130: CrystalQuick Plus: square well design, flat bottom, 3 wells per reservoir, hydrophobic surface</p> <p>609171: CrystalQuick: square well design, flat bottom, 1 well per reservoir, low profile</p> <p>609180: CrystalQuick Plus: square well design, flat bottom, 1 well per reservoir, low profile, hydrophobic surface</p> <p>Appropriate lids for hanging drop (CrystalDrop), Cat. No. 609050 / 609150</p>
1.2	Dimensions	See Customer Drawing
1.3	Volume per well	<p>609101, -130: max. volume per well: 4.1 µl max. volume per reservoir: 320 µl</p> <p>609120: max. volume per well: 1.9 µl max. volume per reservoir: 320 µl</p> <p>609171, -180: max. volume per well: 3.9 µl max. volume per reservoir: 140 µl</p>
1.4	Material / Resin	PS (Polystyrene), free of heavy metal
1.5	Colour	Clear
1.6	Sterilization	No
1.7	Quality Control	<p>- <u>Raw Material-Control</u>: physical testing</p> <p>- <u>Product-Control</u>: testing of attributive and variable characteristics in accordance with the valid specification</p>
1.8	Intended Use	General laboratory product for the processing and storage of samples to be used by qualified personnel in a laboratory environment.
1.9	Other Information	For single use only

2.	Features	
2.1	Basic features	609130, -180: hydrophobic surface
2.2	Temperature range	For application: -20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	609101, -120, -130	609171, -180
3.1	Pieces / Bag	10	20
3.2	Pieces / Box	40	80
3.3	Lot-No.	E JJ MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	E JJ MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-	-

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 04	Date 3 February 2023	Date 13 April 2023	Date 13 April 2023	
Date 02.12.2014	Name S. Kaelberer	Name P. Wachter	Name Dr. Chen-Ket Chai	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.