


PDS No. 6281xx / 628160-TRI	<b>PRODUCT DATA SHEET</b>				Page 1 of 1
Revision 08	Petri Dish, 60 x 15 mm				
	Greiner Item-No. 6281xx / 628160-TRI				
Valid for Item-No.:	<b>628102</b>	<b>628160</b> (sterile)	<b>628161</b> (sterile)	<b>628160-TRI</b> (sterile)	

1.	Description / Specification	
1.1	Description	Petri Dish with vents 628102: standard 628160: CELLSTAR®, physical surface treatment, sterile 628161: sterile 628160-TRI: CELLSTAR®, physical surface treatment, sterile, triple packed
1.2	Dimensions	See Customer Drawing Total weight: 7.4 – 7.8 g
1.3	Volume	Max. volume: 28 ml Working volume: 6-7 ml 628160, -160-TRI: growth area: 21 cm <sup>2</sup>
1.4	Material / Resin	Dish: PS (Polystyrene), free of heavy metal Lid: PS (Polystyrene), free of heavy metal
1.5	Colour	Dish: clear Lid: clear
1.6	Sterilisation	628102: no 628160, -161: SAL 10 <sup>-3</sup> 628160-TRI: SAL 10 <sup>-6</sup>
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	628102, -161: general laboratory product for bacteriology to be used by qualified personnel in a laboratory environment. 628160, -160-TRI: general laboratory product for cell culture to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- For single use only - Elevations for ventilation of culture

2.	Features	
2.1	Basic features	628102, 161: - 628160, -160-TRI: free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	628102: n/a 628160, -160-TRI: 4 years 628161: 5 years
2.7	Other Information	-

3.	Packaging	628102, -161	628160	628160-TRI
3.1	Pieces / Bag	20	10	10 (triple packed)
3.2	Pieces / Box	600	600	300
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)		
3.4	Other Information	628160, -160-TRI: Certificate of Quality to download		

4.	Other Information
4.1	628160 / 628160-TRI: Research use only. Not for diagnostics.

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 07	Date 5 September 2022	Date 14 September 2022	Date 15 September 2022	
Date 26.08.2015	Name S. Kaelberer	Name Dr. R. Daum	Name A. Illig	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.