

|                         |                                                |           |        |           |                                                                                     |  |
|-------------------------|------------------------------------------------|-----------|--------|-----------|-------------------------------------------------------------------------------------|--|
| PDS No. 78407x / 784101 | PRODUCT DATA SHEET                             |           |        |           | Page 1 of 1                                                                         |  |
| Revision 07             | 384 Well Microplate, PS, Small Volume™, HiBase |           |        |           |  |  |
|                         | Item-No. 78407x / 784101                       |           |        |           |                                                                                     |  |
| Valid for Item-No.:     | 784075                                         | 784075-25 | 784076 | 784076-25 | 784101                                                                              |  |

| 1.  | Description / Specification |                                                                                                                                                           |
|-----|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1.1 | Description                 | PS Microplate, 384 well, solid bottom, Small Volume™, HiBase<br>784101: no treatment<br>784075, 784075-25, 784076, 784076-25: medium binding              |
| 1.2 | Dimensions                  | See Customer Drawing                                                                                                                                      |
| 1.3 | Volume per well             | Total volume: 28 µl (mathematical calculated)<br>Working volume: 4 – 25 µl                                                                                |
| 1.4 | Material / Resin            | PS (Polystyrene), free of heavy metal                                                                                                                     |
| 1.5 | Colour                      | 784075, 784075-25: white<br>784076, 784076-25: black<br>784101: clear                                                                                     |
| 1.6 | Sterilization               | No                                                                                                                                                        |
| 1.7 | Quality Control             | Raw Material-Control: physical testing<br>Product-Control: testing of attributive and variable characteristics in accordance with the valid specification |
| 1.8 | Other Information           | For single use only                                                                                                                                       |

| 2.  | Features                 |                                                                                                                                                                    |
|-----|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 2.1 | Basic features           | Free of detectable DNase/RNase, human DNA and pyrogens                                                                                                             |
| 2.2 | Temperature range        | -20°C to +60°C                                                                                                                                                     |
| 2.3 | Autoclavability          | No                                                                                                                                                                 |
| 2.4 | Centrifugation, max. RCF | 800 x g: swinging-bucket rotor                                                                                                                                     |
| 2.5 | Chemical Resistance      | See homepage:<br><a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a> |
| 2.6 | Shelf life               | 784075, 784075-25, 784076, 784076-25: 4 years<br>784101: n/a                                                                                                       |
| 2.7 | Other Information        | -                                                                                                                                                                  |

| 3.  | Packaging         | 784075, -076, -101                                                     | 784075-25 / 784076-25                                                  |
|-----|-------------------|------------------------------------------------------------------------|------------------------------------------------------------------------|
| 3.1 | Pieces / Bag      | 10                                                                     | 25                                                                     |
| 3.2 | Pieces / Box      | 40                                                                     | 150                                                                    |
| 3.3 | Lot-No.           | E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.) | E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.) |
| 3.4 | Other Information | Certificate of Quality                                                 | Certificate of Quality                                                 |

| 4. | Other Information |
|----|-------------------|
|    | -                 |

Data Sheet subject to change without notice!

| Prior Issue        | Drawn                 | Approved              | Released              | <b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved. |
|--------------------|-----------------------|-----------------------|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Revision 06        | Date<br>14 March 2025 | Date<br>14 March 2025 | Date<br>14 March 2025 |                                                                                                                                                                                                                                                                                                       |
| Date<br>14.09.2020 | Name<br>P.Jünger      | Name<br>P.Wachter     | Name<br>A.Mackowski   |                                                                                                                                                                                                                                                                                                       |
|                    |                       |                       |                       |                                                                                                                                                                                                                                                                                                       |

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.