


PDS No. 95x710	<b>PRODUCT DATA SHEET</b>				Page 1 of 1
Revision 00	Media Bottle with Standard Screw Cap				
	Greiner Item-No. 95x710				
Valid for Item-No.:	950710 (sterile)	951710 (sterile)	952710 (sterile)	953710 (sterile)	

1.	Description / Specification	
1.1	Description	Squared media bottle (tetragonal design) with standard screw cap, scale engraving, sterile 950710: 500 ml 951710: 100 ml 952710: 1000 ml 953710: 250 ml
1.2	Dimensions	See Customer Drawing
1.3	Volume	950700: 500 ml 951700: 100 ml 952700: 1000 ml 953710: 250 ml
1.4	Material / Resin	<u>Bottle</u> : PET (Polyethylene Terephthalate), free of heavy metal <u>Screw cap</u> : HDPE (High Density Polyethylene), free of heavy metal
1.5	Colour	<u>Bottle</u> : transparent <u>Screw cap</u> : translucent white
1.6	Sterilization	SAL 10 <sup>-6</sup>
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for the storage of media and buffers to be used by qualified personnel in a laboratory environment.
1.9	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -20°C to +40°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	5 years
2.7	Other Information	-

3.	Packaging	950710	951710	952710	953710
3.1	Pieces / bag	25	49	12	36
3.2	Pieces / Box	50	98	24	108
3.3	Lot-No.	7-digit SAP number	7-digit SAP number	7-digit SAP number	7-digit SAP number
3.4	Other Information	-	-	-	-

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision	Date	Date	Date	
-	6 March 2023	12 June 2023	12 June 2023	
Date	Name	Name	Name	
-	S. Kaelberer	Dr. T. Schreiber	Dr. Chen-Ket Chai	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.