


PDS No. 65020x	<b>PRODUCT DATA SHEET</b>				Page 1 of 1
Revision 03	96 Well Microplate, PP, U-bottom, Colored				
	Greiner Item-No. 65020x				
Valid for Item-No.:	<b>650203</b>	<b>650204</b>	<b>650205</b>	<b>650206</b>	

1.	Description / Specification	
1.1	Description	PP Microplate, 96 well, solid U-bottom, chimney well, alphanumeric well coding, different colours
1.2	Dimensions	See Customer Drawing
1.3	Volume per well	Total volume: 355 µl (mathematical calculated) Working volume: 50 - 300 µl
1.4	Material / Resin	PP (Polypropylene), free of heavy metal
1.5	Color	650203: red 650204: blue 650205: green 650206: yellow
1.6	Sterilization	No
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for the processing and storage of samples to be used by qualified personnel in a laboratory environment.
1.9	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens.
2.2	Temperature range	For application: -196°C to +121°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	10
3.2	Pieces / Box	100
3.3	Lot-No.	B YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality to download

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 02	Date 8 February 2022	Date 27 May 2022	Date 30 May 2022	
Date 27.11.2014	Name S. Kaelberer	Name P. Wachter	Name A. Illig	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.