


PDS No. 6626xx	PRODUCT DATA SHEET				Page 1 of 1
Revision 05	ThinCert® Cell Culture Inserts for 24 Well Plates				
	Greiner Item-No. 6626xx				
Valid for Item-No.:	662610 (sterile)	662630 (sterile)	662631 (sterile)	662638 (sterile)	662640 (sterile)
	662641 (sterile)				

1.	Description / Specification																						
1.1	Description	ThinCert® Cell Culture Inserts, hanging geometry, sealed capillary pore membrane, physical surface treatment, sterile																					
		<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 33%;">Pore size (µm)</th> <th style="width: 33%;">Membrane thickness (µm)</th> <th style="width: 33%;">Pore density (cm²)</th> </tr> </thead> <tbody> <tr> <td>662610: 1,0</td> <td>22 +/- 3</td> <td>662610: 2 x 10⁶</td> </tr> <tr> <td>662630: 3,0</td> <td>20 +/- 3</td> <td>662630: 0,6 x 10⁶</td> </tr> <tr> <td>662631: 3,0</td> <td>20 +/- 3</td> <td>662631: 2 x 10⁶</td> </tr> <tr> <td>662638: 8,0</td> <td>15 +/- 3</td> <td>662638: 0,15 x 10⁶</td> </tr> <tr> <td>662640: 0,4</td> <td>22 +/- 3</td> <td>662640: 1 x 10⁸</td> </tr> <tr> <td>662641: 0,4</td> <td>22 +/- 3</td> <td>662641: 2 x 10⁶</td> </tr> </tbody> </table>	Pore size (µm)	Membrane thickness (µm)	Pore density (cm ²)	662610: 1,0	22 +/- 3	662610: 2 x 10 ⁶	662630: 3,0	20 +/- 3	662630: 0,6 x 10 ⁶	662631: 3,0	20 +/- 3	662631: 2 x 10 ⁶	662638: 8,0	15 +/- 3	662638: 0,15 x 10 ⁶	662640: 0,4	22 +/- 3	662640: 1 x 10 ⁸	662641: 0,4	22 +/- 3	662641: 2 x 10 ⁶
Pore size (µm)	Membrane thickness (µm)	Pore density (cm ²)																					
662610: 1,0	22 +/- 3	662610: 2 x 10 ⁶																					
662630: 3,0	20 +/- 3	662630: 0,6 x 10 ⁶																					
662631: 3,0	20 +/- 3	662631: 2 x 10 ⁶																					
662638: 8,0	15 +/- 3	662638: 0,15 x 10 ⁶																					
662640: 0,4	22 +/- 3	662640: 1 x 10 ⁸																					
662641: 0,4	22 +/- 3	662641: 2 x 10 ⁶																					
1.2	Dimensions	See customer drawing Culture surface (growth area): 33,6 mm ²																					
1.3	Volume	Working volume ThinCert®: 0,1 – 0,35 ml Working volume well: 0,4 – 1,2 ml																					
1.4	Material / Resin	ThinCert® inserts: PS (Polystyrene), free of heavy metal Capillary pore membrane: PET, free of heavy metal																					
1.5	Colour	ThinCert® inserts: clear Capillary pore membrane: 662610, -630, -641: transparent 662631, -638, -640: translucent																					
1.6	Sterilisation	SAL 10 ⁻³																					
1.7	Quality Control	<u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification																					
1.8	Intended Use	General laboratory product for cell culture to be used by qualified personnel in a laboratory environment																					
1.9	Other Information	USP class VI certified																					

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	4 years
2.7	Other Information	-

3.	Packaging	ThinCert® inserts	Multiwell plates
3.1	Pieces (single packed; paper/plastic)	1	1
3.2	Pieces / Box	48	2
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	
3.4	Other Information	Certificate of Quality to download	

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 04	Date 9 August 2023	Date 10 August 2023	Date 10 August 2023	
Date 23.06.2021	Name S. Kaelberer	Name Dr. T. Binder	Name Dr. C.-K. Chai	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.