


PDS No. 6701xx	<b>PRODUCT DATA SHEET</b>			Page 1 of 1
Revision 05	<b>CELLSTAR® OneWell Plate</b>			
	Greiner Item-No. 6701xx			
Valid for Item-No.:	<b>670102</b>	<b>670180</b> (sterile)	<b>670190</b> (sterile)	

1.	Description / Specification	
1.1	Description	Cell culture plate with 1 rectangular chamber, vents and lid (single position, high profile) 670102: non-surface treatment, non-sterile 670180: physical surface treatment, sterile 670190: non-surface treatment, sterile
1.2	Dimensions	See Customer Drawing
1.3	Volume	Total volume: 113.7 ml 670180: growth area: 95.0 cm <sup>2</sup>
1.4	Material / Resin	Plate and lid: PS (Polystyrene), free of heavy metal
1.5	Colour	Plate and lid: clear
1.6	Sterilisation	670102: - 670180, -190: SAL 10 <sup>-3</sup>
1.7	Quality Control	<u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for cell culture to be used by qualified personnel in a laboratory environment.
1.9	Other Information	For single use only

2.	Features	
2.1	Basic features	670102: - 670180, -190: free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	670102: n/a 670180: 4 years 670190: 5 years
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	8
3.2	Pieces / Box	32
3.3	Lot-No.	M JJ ZZZZ (manufacturing facility, year, consecutive -No.)
3.4	Other Information	-

4.	Other Information
4.1	Research use only. Not for diagnostics.

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 04	Date 4 April 2022	Date 4 October 2022	Date 5 October 2022	
Date 22.04.2015	Name S. Kaelberer	Name P. Wachter	Name A. Illig	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.