


PDS No. 662145	PRODUCT DATA SHEET	Page 1 of 1
Revision 04	CrystalBridge™	 greiner bio-one
	Greiner Item-No. 662145	

1.	Description / Specification	
1.1	Description	CrystalBridge™ insert, concave well bottom for sitting drop experiments, fitting exactly into the wells of the ComboPlates™ (Cat. No. 662050 / 662150 / 662850)
1.2	Dimensions	See customer drawing
1.3	Volume	Max. volume: 45 µl
1.4	Material / Resin	PS (Polystyrene), free of heavy metal
1.5	Colour	Clear
1.6	Sterilization	No
1.7	Quality Control	<u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	-
2.2	Temperature range	-20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	250
3.2	Pieces / Box	Quantity on request
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-

4.	Other Information	
	-	

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 03	Date 12 October 2018	Date 12 October 2018	Date 12 October 2018	
Date 07.03.2018	Name S. Kaelberer	Name Dr. S. Sigle	Name A. Illig	

DISCLAIMER: The description of a certain product can only be considered as a guidance because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.