

PRODUCT DATA SHEET

Petri Dish 90x15mm

Greiner item:	630164	630165
---------------	--------	--------

1.	Description/Specification	
1.1	Description	630164: 90x15mm, sterile 630165: 90x15mm, non sterile
1.2	Dimensions	90x15mm
1.3	Material	Dish: PS (Polystyrene), free of heavy metal Lid: PS (Polystyrene), free of heavy metal
1.4	Colour	Clear
1.5	Sterilisation	630164: SAL 10 ⁻⁶ 630165: Non sterile
1.6	Other information	For single use only

2.	Features	
2.1	Autoclavability	No
2.2	Centrifugatio, max RCF	N/A
2.3	Shelf Life	630164: 05 years after production date. 630165: N/A
2.6	Other information	-

3.	Packaging	
3.1	Pieces / Bag	20
3.2	Pieces / Box	960
3.3	Lot-No	F YY MM 00A (manufacturing facility, year, month, consecutive N°).
3.4	Box dimensions (cm)	59 x 76 x 33cm
3.5	Other information	-