


PDS No. 661960/ 661975	<b>PRODUCT DATA SHEET</b>		Page 1 of 1
Revision 04	Cell Culture Flask, 650 ml, Advanced TC		
	Greiner Item-No. 661960 / 661975		
Valid for Item-No.:	<b>661960</b> (sterile)	<b>661975</b> (sterile)	

1.	Description / Specification	
1.1	Description	Cell Culture Flask, 650 ml, high flask design, canted neck, printed graduation (scale 25-200 ml), writing area, sterile <i>661960</i> : Standard Cell Culture Flasks, Advanced TC surface <i>661975</i> : Filter Cap Cell Culture Flasks, Advanced TC surface
1.2	Dimensions	<u>Flask</u> : see Customer Drawing <i>661975</i> : pore size of filter membrane: 0.2 µm
1.3	Volume	Total volume: 650 ml Working volume: 20-85 ml Growth area: 175 cm <sup>2</sup>
1.4	Material / Resin	<u>Flask</u> : PS (Polystyrene), free of heavy metal <u>Cap</u> : HDPE (High Density Polyethylene), free of heavy metal <u>Filter</u> : PET (Polyethylene Terephthalate)
1.5	Colour	<u>Flask</u> : clear; print: graduation and writing area white <u>Cap</u> : blue <u>Filter</u> : white
1.6	Sterilization	SAL 10 <sup>-3</sup>
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for cell culture to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- For single use only - Expiry date and Lot-No. printed on bottom of flask

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: +4°C to +37 °C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	2 years
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	4
3.2	Pieces / Box	40
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality to download

4.	Other Information
4.1	Research use only. Not for diagnostics.

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 03	Date 17 December 2020	Date 18 January 2021	Date 25 February 2021	
Date 22.04.2015	Name S. Kaelberer	Name R. Daum	Name A. Illig	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.