

|                     |                                                                   |                  |  |  |                                                                                     |  |
|---------------------|-------------------------------------------------------------------|------------------|--|--|-------------------------------------------------------------------------------------|--|
| PDS No. 6278xx      | PRODUCT DATA SHEET                                                |                  |  |  | Page 1 of 1                                                                         |  |
| Revision 07         | CELLview, Cell Culture Dish, 35 x 10 mm,<br>Glass Bottom, Sterile |                  |  |  |  |  |
|                     | Greiner Item-No. 6278xx                                           |                  |  |  |                                                                                     |  |
| Valid for Item-No.: | 627860 (sterile)                                                  | 627870 (sterile) |  |  |                                                                                     |  |

| 1.  | Description / Specification |                                                                                                                                                                                                                                                                                                                         |
|-----|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1.1 | Description                 | Cell culture dish with glass bottom for microscopic applications, with vents, sterile<br><u>Glass bottom</u> : maximal spectral transmission, no autofluorescence, no depolarisation of light<br>627860: with one compartment, physical surface treatment<br>627870: with four compartments, physical surface treatment |
| 1.2 | Dimensions                  | See Customer Drawing<br><u>Weight</u> : 627860 lid: 1.2 – 1.5 g / dish: 1.6 – 1.8 g<br>627870 lid: 1.2 – 1.5 g / dish: 2.5 – 2.8 g<br>Glass thickness: 0.17 mm (+/- 0.01 mm)                                                                                                                                            |
| 1.3 | Volume                      | 627860: working volume: 5.0 ml<br>627870: working volume: 0.1-0.5 ml / per compartment<br>627860: growth area: 8.7 cm <sup>2</sup><br>627870: growth area: 1.9 cm <sup>2</sup> / per compartment                                                                                                                        |
| 1.4 | Material / Resin            | <u>Dish</u> : PS (Polystyrene), free of heavy metal<br><u>Bottom</u> : borosilicate glass; hydrolytic class 1 (DIN-ISO 719)<br><u>Lid</u> : PS (Polystyrene), free of heavy metal                                                                                                                                       |
| 1.5 | Colour                      | <u>Dish</u> : clear<br><u>Glass Bottom</u> : high transparent<br><u>Lid</u> : clear                                                                                                                                                                                                                                     |
| 1.6 | Sterilisation               | SAL 10 <sup>-3</sup>                                                                                                                                                                                                                                                                                                    |
| 1.7 | Quality Control             | - <u>Raw Material-Control</u> : physical testing<br>- <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification                                                                                                                                           |
| 1.8 | Intended Use                | General laboratory product for cell culture to be used by qualified personnel in a laboratory environment                                                                                                                                                                                                               |
| 1.9 | Other Information           | - For single use only<br>- Manufactured according to ISO 8255-1:2011                                                                                                                                                                                                                                                    |

| 2.  | Features                 |                                                                                                                                                                              |
|-----|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 2.1 | Basic features           | Free of detectable DNase/RNase, human DNA and pyrogens.<br>Contents non-cytotoxic<br><u>Adhesive</u> : resistant to DMSO and 70% ethanol, compatible with cell culture media |
| 2.2 | Temperature range        | For application: + 4°C to +37°C                                                                                                                                              |
| 2.3 | Autoclavability          | No                                                                                                                                                                           |
| 2.4 | Centrifugation, max. RCF | N/A                                                                                                                                                                          |
| 2.5 | Chemical Resistance      | See homepage:<br><a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>           |
| 2.6 | Shelf life               | 2 years                                                                                                                                                                      |
| 2.7 | Other Information        | -                                                                                                                                                                            |

| 3.  | Packaging         |                                                                        |
|-----|-------------------|------------------------------------------------------------------------|
| 3.1 | Pieces / Bag      | 10                                                                     |
| 3.2 | Pieces / Box      | 40                                                                     |
| 3.3 | Lot-No.           | E JJ MM XXX (manufacturing facility, year, month, consecutive SAP-No.) |
| 3.4 | Other Information | Certificate of Quality to download                                     |

| 4.  | Other Information                       |
|-----|-----------------------------------------|
| 4.1 | Research use only. Not for diagnostics. |

Data Sheet subject to change without notice!

| Prior Issue     | Drawn              | Approved           | Released           | <b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved. |
|-----------------|--------------------|--------------------|--------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Revision 06     | Date 14 March 2025 | Date 14 March 2025 | Date 14 March 2025 |                                                                                                                                                                                                                                                                                                       |
| Date 15.09.2022 | Name P.Jünger      | Name P.Wachter     | Name A.Mackowski   |                                                                                                                                                                                                                                                                                                       |

**DISCLAIMER:** The description of a certain product can only be considered as a guidance because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.